


JOBST® Elvarex®

FASE 2 CDP


Su misura & Standard

Modelli **Su misura** tutori elastocontenitivi **arti inferiori**




CCL1 18-21 mmHg	•	•	•	•	•	•	•	•	•
CCL2 25-32 mmHg	•	•	•	•	•	•	•	•	•
CCL3 circa 40 mmHg	•	•	•	•	•	•	•		
CCL3 forte circa 44 mmHg	•	•	•	•	•	•	•		
CCL4 49-70 mmHg	•	•	•						
CCL4 Super 60-90 mmHg	•	•	•						

Varianti




Modelli **Standard** tutori elastocontenitivi **arti inferiori**



CCL1 18-21 mmHg	•	•
CCL2 25-32 mmHg	•	•

Tabella delle misure

PUNTO DI MISURAZIONE		4	5-6	7-8	9-10	11-12	13-14
Misure <i>c</i> Circonferenza cm	<i>cG</i>	42 - 44	46 - 48	50 - 52	54 - 56	58 - 60	62 - 64
	<i>cF</i>	39 - 41	42 - 44	45 - 47	48 - 50	51 - 53	54 - 56
	<i>cE</i>	29 - 31	32 - 34	35 - 37	38 - 40	41 - 43	44 - 46
	<i>cD</i>	26 - 28	29 - 31	32 - 34	35 - 37	38 - 40	41 - 43
	<i>cC</i>	29 - 30	31 - 33	34 - 36	37 - 39	40 - 42	43 - 45
	<i>cB1</i>	23 - 24	25 - 26	27 - 28	29 - 30	31 - 32	33 - 34
	<i>cY</i>	27 - 28	29 - 30	31 - 32	33 - 34	35 - 36	37 - 38
	<i>cA/cB</i>	18 - 19	20 - 21	22 - 23	24 - 25	26 - 27	28 - 29



AD: /D normale = 40-43 cm, /D corto = 36-39 cm • **AG:** /G = 69-74 cm • **AG-T:** /K = 70-75 cm